

**General Purpose Homogenisers, Borosilicate glass****GPE****General Purpose Homogenisers, Borosilicate glass**

- Precision Bore Borosilicate glass mortars with a range from 5.0 to 55ml capacity
- Borosilicate glass pestle
- All parts are fully autoclavable

Grinding Clearance Pestle and Mortar: 0.152 to 0.254mm



Code	Description	Packaging
LLG06264841	Homogenisers 55 ml, Chamber Length 125 mm, Mortar Overall length 175 mm, Pestle Length 295 mm, Pestle Diam. 25 mm	1 pz.
LLG06264838	Homogenisers 10 ml, Chamber Length 75 mm, Mortar Overall length 115 mm, Pestle Length 215 mm, Pestle Diam. 13 mm	1 pz.
LLG06264839	Homogenisers 15 ml, Chamber Length 80 mm, Mortar Overall length 120 mm, Pestle Length 220 mm, Pestle Diam. 15 mm	1 pz.
LLG06264840	Homogenisers 30 ml, Chamber Length 105 mm, Mortar Overall length 150 mm, Pestle Length 250 mm, Pestle Diam. 19 mm	1 pz.
LLG06264837	Homogenisers 5 ml, Chamber Length 65 mm, Mortar Overall length 95 mm, Pestle Length 195 mm, Pestle Diam. 10 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

