

**Accessories for Kjeldahl rapid sample digestors****behr****Accessories for Kjeldahl rapid sample digestors**

behrotest® standard digestion vessels have standard dimensions. That means they fit into digestion and distillation units of most brands. Of course, it also means that standard digestion vessels of most brands can be used with behr distillation units.



Code	Description	Packaging
LLG09843753	Accessories for Kjeldahl rapid sample digestors, Type AFS, Antifoam tablets	1000 pz.
LLG09843969	Accessories for Kjeldahl rapid sample digestors, Type KT 4, Catalyst tablets (3.5 g K ₂ SO ₄ ; 3.5 mg Se)	1 pz.
LLG09843754	Accessories for Kjeldahl rapid sample digestors, Type WP, Nitrogen-free weighing boats	250 pz.
LLG09843790	Accessories for Kjeldahl rapid sample digestors, Type EG 4/ 500, Vessel holding frame for four 500ml vessels. Component of InKjel 450 M and P	
LLG09843766	Accessories for Kjeldahl rapid sample digestors, Type SR 3i, Sample vessel, 250ml, for InKjel*	1 pz.
LLG09843775	Accessories for Kjeldahl rapid sample digestors, Type KJ 500, Sample vessel, 500ml, for InKjel*	1 pz.
LLG09843751	Accessories for Kjeldahl rapid sample digestors, Type SR 5, Sample vessel 500ml, for K8 B*	1 pz.
LLG09843776	Accessories for Kjeldahl rapid sample digestors, Type KJ 750, Sample vessel, 750ml, for InKjel*	1 pz.
LLG09843767	Type KT 1, Catalyst tablets (5.0 g K ₂ SO ₄ ; 0.5 g CuSO ₄)	1000 pz.
LLG09843752	Accessories for Kjeldahl rapid sample digestors, Type KT 3, Catalyst tablets (3.5 g K ₂ SO ₄ ; 0.4 g CuSO ₄)	1000 pz.
LLG09843791	Accessories for Kjeldahl rapid sample digestors, Type EG 4/ 750, Vessel holding frame for four 750ml vessels. Component of InKjel 475 M and P	
LLG09843788	Accessories for Kjeldahl rapid sample digestors, Type EG 6, Vessel holding frame for six 250ml vessels. Component of InKjel 625 M and P	
LLG09843789	Accessories for Kjeldahl rapid sample digestors, Type EG 12, Vessel holding frame for twelve 250ml vessels. Component of InKjel 1225 M and P	
LLG09843773	Accessories for Kjeldahl rapid sample digestors, Type KT 2, Catalyst tablets (5.0 g K ₂ SO ₄ ; 0.15 g CuSO ₄ ; 0.15 g TiO ₂)	1000 pz.
LLG09920710	Boiling stones behrotest®, 100 g	1 pz.
LLG09843750	Accessories for Kjeldahl rapid sample digestors, Type SR 4, Sample vessel 100ml, for InKjel*	1 pz.



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

