

**Magnetic stirrer MIXdrive 1 XL with external control****43953****Magnetic stirrer MIXdrive 1 XL with external control**

Due to the external control, these magnetic stirrers are suitable for use in ovens, incubators, safety cabinets or for underwater operation, depending on the device. Wide speed ranges make them ideal for viscous media. Low space requirement due to the flat and compact design with maximum stirring surface.

- Maintenance-free 2mag motor drive
- Fully sealed housing
- Extremely robust construction and design
- Easy to clean and maximum resistance
- 3 years warranty
- Made in Germany

**Please always order control units separately.**

**Magnetic stirrer for large quantities of stirring and with a very powerful magnetic coupling for chemical and biotechnical applications.**

**-High power transmission over longer distances or for thick vessel bottoms**

**-Optimal for use in ovens up to 100 °C**

**-Stirring bar monitoring mixWATCH in combination with MIXcontrol XL**

**Specifications**

**Stirring quantity (H<sub>2</sub>O): 1 ... 40 l**

**Permissible ambient temperature: 5 ... 100 °C at 95 % r. H.**

**Housing: stainless steel**

**Dimesnions (W x D x H): 225 x 225 x 38 mm**

**Weight: approx. 3.4 kg**

**IP code: IP 66**



Code	Description	Packaging
LLG09646029	Control unit MIXcontrol 1XL 60-2000 rpm, 20 W, IP20, 200x155x48mm	1 pz.
LLG09646028	Magnetic stirring system MIXdrive 1XL (+150°C), 60-2000rpm, IP66, external control, 48 VDC	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

