



## Complete compact extraction systems, with heating

### Complete compact extraction systems, with heating

Complete compact extraction systems with stand, heating module, mount, hoses and glassware (round-bottom reaction flask, extractor, Dimroth condenser for extraction). Infinitely variable heating regulation.

- practical stand for holding the condensers securely between extractions
- extractor sizes from 30 ml to 5000 ml
- the use of extractors with a stopcock eliminates the need for any additional distillation after extraction
- Condensers with screwed-on tubing nipples



Code	Description	Packaging
LLG09843778	Complete compact extraction systems, with heating, Type KEX 100, for 100ml extraction with 250 ml distilling flasks	1 pz.
LLG09843786	Complete compact extraction systems, with heating, Type KEX 100 F, for 100ml extraction with stopcock and 250 ml distilling flasks	1 pz.
LLG09843779	Complete compact extraction systems, with heating, Type KEX 250 F, for 250ml extraction with stopcock and 500 ml distilling flasks	1 pz.
LLG09843777	Complete compact extraction systems, with heating, Type KEX 30, for 30ml extraction with 100 ml distilling flasks	1 pz.
LLG09843836	Complete compact extraction systems, with heating, Type KEX 30 F, for 30ml extraction with stopcock and 100 ml distilling flasks	1 pz.
LLG09843835	Complete compact extraction systems, with heating, Type KEX 60, for 60ml extraction with 250 ml distilling flasks	1 pz.
LLG09843837	Complete compact extraction systems, with heating, Type KEX 60 F, for 60ml extraction with stopcock and 250 ml distilling flasks	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

