

**Aspirator bottles, series 350, HDPE****Kautex****Aspirator bottles, series 350, HDPE**

Food-grade product according to Regulation (EC) No 10/2011. Natural with printed blue graduations (in L) and screwed closure. With or without threaded connection at base. Inside neck diameter 50mm.

Please order stopcock separately.



Code	Description	Packaging
LLG09011010	Aspirator bottles, series 350, HDPE, 10 Litres	1 pz.
LLG09011071	Aspirator bottles, series 350, HDPE, 10 Litres	1 pz.
LLG09011025	Aspirator bottles, series 350, HDPE, 25 Litres	1 pz.
LLG09011072	Aspirator bottles, series 350, HDPE, 25 Litres	1 pz.
LLG09011009	Aspirator bottles, series 350, HDPE, 5 Litres	1 pz.
LLG09011069	Aspirator bottles, series 350, HDPE, 5 Litres	1 pz.
LLG09011060	Aspirator bottles, series 350, HDPE, 60 Litres	1 pz.
LLG09011076	Aspirator bottles, series 350, HDPE, 60 Litres	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

