



Storage and transport containers, PP / HDPE

Storage and transport containers, PP / HDPE

White. Made of thermoplastic resin. Stackable.

The raw material used is physiologically compatible and corresponds with the directives of the Federal Board of Health.



Code	Description	Packaging
LLG09301641	Storage and transport containers, PP / HDPE, 10 Litres	1 pz
LLG09301663	Storage and transport containers, PP / HDPE, 23 Litres	1 pz
LLG09301647	Storage and transport containers, PP / HDPE, 28 Litres	1 pz
LLG09301645	Storage and transport containers, PP / HDPE, 31 Litres	1 pz
LLG09301661	Storage and transport containers, PP / HDPE, 35 Litres	1 pz
LLG09301649	Storage and transport containers, PP / HDPE, 50 Litres	1 pz
LLG09301643	Storage and transport containers, PP / HDPE, 9 Litres	1 pz
LLG09301651	Storage and transport containers, PP / HDPE, 96 Litres	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

