

**Dry Milk Butyrometer Teichert with Accessories****Gerber****Dry Milk Butyrometer Teichert with Accessories**

Weighing method for 2,5 g dry milk, closed bulb with frosted label, with rubber stopper and glass nail.



Code	Description	Packaging
LLG09112569	Dry Milk Butyrometer Teichert with Accessories, Fat 0 - 35 %, Graduation 0.5 ( $\pm 0.25$ )	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

