



Accessories for Butyrometers

Accessories for Butyrometers



Code	Description	Packaging
LLG09112573	Accessories for Butyrometers, GERBAL-M stopper with bush, length 23.5mm, Accessories Code A	1 pz.
LLG09112578	Glass Beaker with 16 holes for Cheese Butyrometers	1 pz.
LLG06262314	Glass beaker with two holes for butter butyrometers	1 pz.
LLG09112577	Glass beaker without holes for Cream Butyrometers	1 pz.
LLG06234037	Glassnail for powder milk butyrometer	1 pz.
LLG09112574	Accessories for Butyrometers, Regulating pin for GERBAL-M stopper, nickel plated brass, Accessories Code B	1 pz.
LLG06202638	Accessories for Butyrometers, Single conical rubber stopper 35x11/16mm, acid-proof for milk, skim milk and Koehler-Cream-Butyrometer, Accessories Code F	1 pz
LLG09112575	Single conical rubber stopper, large, 30 x 17.5/22mm, with hole for beaker	1 pz.
LLG06202636	Single conical rubber stopper, large, 30 x 17.5/22mm, without hole	
LLG09112576	Single conical rubber stopper, small, 15x10.5/13 mm	1 pz.
LLG06264750	Weighing funnel of glass, for weighing and filling the milkpowder into the butyrometer	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

