

**Water jet pumps, PP****BRAND****Water jet pumps, PP**

With constant ultimate vacuum, high suction flow rate and very low water consumption. Can be connected to mains water system in a number of different ways using the adapter supplied and reducing adapters that are available as accessories. For continuous use at temperatures up to 80°C max. High chemical resistance as the media being pumped only comes into contact with PP, FKM and PTFE. Integral non-return valve increases safety.

Comprises: Water jet pump, including:

Mains water connections: R ¾" sleeve nut, R ½" reducing adapter and flexible tubing connection (nozzle) with external diameter from 10mm to 12mm.

Vacuum connection: Detachable nozzle with external diameter from 6mm to 9mm and GL 14 screw cap.

Performance characteristics at 3.5 bar water supply pressure at 12°C

Water consumption: 190L/hr.

Ultimate vacuum: 16 mbar

Flow rate at atmospheric pressure: 400L/hr. free air



Code	Description	Packaging
LLG07020037	Water jet pumps, PP, Type Reducing adapter R 3/8" for water jet filter pump	1 pz
LLG07020038	Water jet pumps, PP, Type Reducing adapter M 22x1 for water jet filter pump (thread for screen tap)	1 pz
LLG09303125	Water jet pumps, PP, Type Water jet pump	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

