

**Exhaust filters V3.0****S.C.A.T.****Exhaust filters V3.0**

The exhaust air filters are equipped with three different activated carbon layers, which bind acids and alkalis and adsorb organic solvent vapours. Depending on the required capacity, filter with different service lives are available. Filter media with a very high performance spectrum for acids and alkalis Optimised for binding buffer solutions and solvent vapours With splash guard 1500 m<sup>2</sup>/g adsorption surface area



Code	Description	Packaging
LLG04681456	Exhaust filter V3.0, With interchangeable label, Filter dimension M, Lifetime 6 months, Thread 14 GL	
LLG04681457	Exhaust filter V3.0, With change indicator, Filter dimension M, Service life 6 months, Thread 14 GL	
LLG04681458	Exhaust filter V3.0, Spare package, with change indicator, Filter dimension M, Service life 6 months, Thread 14 GL	
LLG04681459	Exhaust filter V3.0, Spare package, with exchange label, Filter dimension M, Service life 6 months, Thread 14 GL	
LLG04681460	Exhaust filter V3.0, With change indicator, Filter dimension S, Service life 3 months, Thread 14 GL	
LLG04681461	Exhaust filter V3.0, Spare package, with change indicator, Filter dimension S, Service life 3 months, Thread 14 GL	
LLG04681462	Exhaust filter V3.0, With change indicator, Filter dimension L, Lifetime 12 months, Thread 14 GL	
LLG04681463	Exhaust filter V3.0, Description Reserve package, with change indicator , Filter dimension L , Service life 12 months , Thread 14 GL	



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

