

**Hand pumps Otal****Bürkle****Hand pumps Otal**

For jerry cans and ballons up to 60l capacity. Universal bungs for container openings of 30 to 70mm diameter.

PP: weak acids, caustic solutions and detergents

PVDF: suitable for acids and alkalis, chlorine bleaching, not for strongly oxidizing acids and caustic solutions

Stainless steel: solvents and liquids of VbF Hazard Class AIII (not easily flammable)



Code	Description	Packaging
LLG09001210	Hand pumps Otal, Diam. 10 mm, Immersion depth 450 mm, Flow rate 4 L / min., Material PP*	1 pz.
LLG09001137	Hand pumps Otal, Diam. 18 mm, Immersion depth 900 mm, Flow rate 26 L / min., Material Stainless steel tubing	1 pz.
LLG09001112	Hand pumps Otal, Diam. 12 mm, Immersion depth 700 mm, Flow rate 14 L / min., Material PP	1 pz.
LLG09001138	Hand pumps Otal, Diam. 12 mm, Immersion depth 700 mm, Flow rate 12 L / min., Material PP / PVC dispensing tube	1 pz.
LLG09001115	Hand pumps Otal, Diam. 15 mm, Immersion depth 700 mm, Flow rate 20 L / min., Material PP	1 pz.
LLG09001135	Hand pumps Otal, Diam. 16 mm, Immersion depth 700 mm, Flow rate 20 L / min., Material PVDF	1 pz.
LLG09001118	Hand pumps Otal, Diam. 18 mm, Immersion depth 700 mm, Flow rate 26 L / min., Material PP	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

