



Beaker tongs, with rubber-sleeved jaws

Bochem

Beaker tongs, with rubber-sleeved jaws

Jaws protected with 2 rubber pads each. Arm length: 160mm.



Code	Description	Packaging
LLG09310010	Beaker tongs, with rubber-sleeved jaws, Length 300 mm, Jaw opening 50 ... 100 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

