

**Microtitre Plate Shakers Titramax 100/101****Heidolph****Microtitre Plate Shakers Titramax 100/101**

-Set and continuously adjust the variable speed on the analog control knob for gentle to powerful mixing

-An analog process timer allows for unattended operation.

Compact and space-saving models for 4 microtiter plates. Titramax 100 with 1.5 mm shaking orbit is ideal for gentle mixing. Titramax 101 with a 3 mm shaking orbit is especially suitable for mixing of viscous media and solids.

Specifications

Max. Load: 2 kg

Platform dimensions (W x D): 220 x 220 mm

Overall dimensions (W x D x H): 245 x 310 x 125 mm



Code	Description	Packaging
LLG09839831	Microtiter plate shaker Titramax 100 max. 4 microtiter plates,orbit 1.5 mm 245x285x125 mm	1 pz.
LLG09839830	Microtiter plates shaker Titramax 1000 complete with 290x258 mm platform,orbit 1.5 mm shaking frequency 100 - 1350 1/min,	1 pz.
LLG09839832	Microtiter plate shaker Titramax 101 max. 4 microtiter plates,orbit 3.0 mm 245x285x125 mm	1 pz.
LLG09839838	Package Titramax -platform shaker Titramax 1000 -heating module Inkubator 1000	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

