

**Centrifuge tubes, round bottom, DURAN****DWK Life Sciences****Centrifuge tubes, round bottom, DURAN®**

To DIN 58970 part 2. Glass Type I / neutral glass as per USP, EP and JP.

Autoclavable.

Maximum loading: 4000xg rcf.



Code	Description	Packaging
LLG09315011	Centrifuge tubes, round bottom, DURAN®, 12 ml, Diam. 16 mm, Height 100 mm	1 pz
LLG09315014	Centrifuge tubes, round bottom, DURAN®, 25 ml, Diam. 24 mm, Height 100 mm	1 pz
LLG09315036	Centrifuge tubes, round bottom, DURAN®, 250* ml, Diam. 56 mm, Height 147 mm	1 pz
LLG09315017	Centrifuge tubes, round bottom, DURAN®, 50 ml, Diam. 34 mm, Height 100 mm	1 pz
LLG09315010	Centrifuge tubes, round bottom, DURAN®, 6 ml, Diam. 12 mm, Height 100 mm	1 pz
LLG09315024	Centrifuge tubes, round bottom, DURAN®, 80* ml, Diam. 40 mm, Height 115 mm	1 pz
LLG09315026	Centrifuge tubes, round bottom, DURAN®, 80 ml, Diam. 44 mm, Height 100 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

