



pH-meter Seven2Go S2

Mettler-Toledo Online

pH-meter Seven2Goä S2

Portable routine model for pH measurements

- New, intuitive menu structure with fast and easy T-Pad navigation in 12 languages
- Robust design with IP67 protection for use in harsh and demanding environments
- Intelligent Sensor Management (ISM®) functionality reduces the risk of errors
- Up to 5 point calibration with automatic buffer recognition and temperature compensation
- 4 predefined buffer groups and 1 user defined buffer group
- Automatic, manual and time dependant endpoint recording
- Data memory storage for 200 data points
- Power requirements: 4 x 1.5V AA or 1.3V NiMH accumulators

Items supplied:

S2-Meter: Instrument, electrode clip, operating instructions, quick guide, test certificate, declaration of conformity and batteries

S2-Standard kit: As S2-Meter, but also with InLab Expert Pro-ISM-IP67, buffer sachets for pH 4.01, 7.00, 9.21, 10.00

S2-Field kit: As S2-Standard kit, but also with bottles for calibration, uGoä carrying case

S2-Food kit: As S2-Field kit, but with InLab Solids Pro-IP67 instead of InLab Expert Pro-ISM-IP67

S2-Light kit: As S2-Standard kit, but with InLab Versatile Pro instead of InLab Expert Pro-ISM-IP67

Specifications**pH**

Measuring range:

Resolution:

Accuracy: -2 ... 20 pH

0.01 / 0.1

±0.01

mV

Measuring range:

Resolution:

Accuracy:



Relative mV-measurement: -1999 ... 1999mV

1

±1

yes

Temperature

Measuring range:

Resolution:

Accuracy: -5.0 ... 105.0 °C

0.1 °C

±0.5 °C

Dimensions (W x D x H): 222 x 70 x 35 mm

Code	Description	Packaging
LLG06265414	pH-meter Seven2Go S2, Type S2-Meter	1 pz
LLG06265415	pH-meter Seven2Go S2, Type S2-Standard kit	1 pz
LLG06265416	pH-meter Seven2Go S2, Type S2-Field kit	1 pz
LLG06265417	pH-meter Seven2Go S2, Type S2-Food kit	1 pz



Code	Description	Packaging
LLG06265418	pH-meter Seven2Go S2, Type S2-Light kit	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

