



Sockets / Cones with NS joint, DURAN tubing

Sockets / Cones with NS joint, DURAN® tubing

Made of DURAN® tubing. DIN 12249 and in addition to DIN standard (*).



Code	Description	Packaging
LLG06222104	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 10/30*, Ext. diam. 8 mm, Length 115 mm	1 pz.
LLG07630287	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 50/42*, Ext. diam. 46 mm, Length 145 mm	1 pz.
LLG07630288	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 55/44*, Ext. diam. 50 mm, Length 145 mm	1 pz.
LLG09011506	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 19/26, Ext. diam. 17 mm, Length 130 mm	1 pz.
LLG06070542	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 5/20*, Ext. diam. 4 mm, Length 115 mm	1 pz.
LLG06203538	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 29/42*, Ext. diam. 26 mm, Length 140 mm	1 pz.
LLG06401432	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 34/45*, Ext. diam. 30 mm, Length 140 mm	1 pz.
LLG06900898	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 14/35*, Ext. diam. 13 mm, Length 140 mm	1 pz.
LLG09011501	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 5/13, Ext. diam. 4 mm, Length 115 mm	
LLG09011502	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 7/16, Ext. diam. 6 mm, Length 115 mm	
LLG09011503	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 10/19, Ext. diam. 8 mm, Length 115 mm	1 pz.
LLG09011504	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 12/21, Ext. diam. 11 mm, Length 130 mm	1 pz.
LLG09011505	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 14/23, Ext. diam. 13 mm, Length 130 mm	1 pz.
LLG09011507	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 24/29, Ext. diam. 22 mm, Length 130 mm	1 pz.
LLG09011508	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 29/32, Ext. diam. 26 mm, Length 130 mm	1 pz.
LLG09011509	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 34/35, Ext. diam. 30 mm, Length 130 mm	
LLG09011510	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 45/40, Ext. diam. 40 mm, Length 145 mm	1 pz.
LLG09011511	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 60/46, Ext. diam. 54 mm, Length 160 mm	
LLG09011512	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 71/51, Ext. diam. 65 mm, Length 180 mm	
LLG09011513	Sockets / Cones with NS joint, DURAN® tubing, Cone, NS 85/55, Ext. diam. 80 mm, Length 180 mm	



Code	Description	Packaging
LLG06072150	Sockets / Cones with NS joint, DURAN® tubing, Cone with drip tip, NS 19/26*, Ext. diam. 17 mm, Length 120 mm	1 pz.
LLG09011545	Sockets / Cones with NS joint, DURAN® tubing, Cone with outlet, NS 14/23, Ext. diam. 13 mm, Length 130 mm	1 pz.
LLG09011546	Sockets / Cones with NS joint, DURAN® tubing, Cone with outlet, NS 19/26, Ext. diam. 17 mm, Length 130 mm	1 pz.
LLG09011547	Sockets / Cones with NS joint, DURAN® tubing, Cone with outlet, NS 24/29, Ext. diam. 22 mm, Length 130 mm	1 pz.
LLG09011548	Sockets / Cones with NS joint, DURAN® tubing, Cone with outlet, NS 29/32, Ext. diam. 26 mm, Length 130 mm	1 pz.
LLG09011550	Sockets / Cones with NS joint, DURAN® tubing, Cone with outlet, NS 45/40, Ext. diam. 40 mm, Length 145 mm	1 pz.
LLG06072920	Sockets / Cones with NS joint, DURAN® tubing, Cone with tip, NS 12/21, Ext. diam. 11 mm, Length 130 mm	1 pz.
LLG06073068	Sockets / Cones with NS joint, DURAN® tubing, Cone with tip, NS 60/46*, Ext. diam. 54 mm, Length 160 mm	1 pz.
LLG06202386	Sockets / Cones with NS joint, DURAN® tubing, Cone with tip, NS 40/38*, Ext. diam. 36 mm, Length 145 mm	1 pz.
LLG06202446	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 40/38*, Ext. diam. 44 mm, Length 130 mm	1 pz.
LLG06303189	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 45/40*, Ext. diam. 48 mm, Length 280 mm	1 pz.
LLG06303190	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 29/32*, Ext. diam. 33 mm, Length 210 mm	1 pz.
LLG09011521	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 5/13, Ext. diam. 8 mm, Length 120 mm	1 pz.
LLG09011522	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 7/16, Ext. diam. 10 mm, Length 120 mm	
LLG09011523	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 10/19, Ext. diam. 13 mm, Length 130 mm	1 pz.
LLG09011524	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 12/21, Ext. diam. 15 mm, Length 130 mm	1 pz.
LLG09011525	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 14/23, Ext. diam. 17 mm, Length 130 mm	1 pz.
LLG09011526	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 19/26, Ext. diam. 22 mm, Length 130 mm	1 pz.



Code	Description	Packaging
LLG09011527	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 24/29, Ext. diam. 28 mm, Length 130 mm	1 pz.
LLG09011528	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 29/32, Ext. diam. 33 mm, Length 130 mm	1 pz.
LLG09011529	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 34/35, Ext. diam. 40 mm, Length 130 mm	
LLG09011530	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 45/40, Ext. diam. 50 mm, Length 145 mm	1 pz.
LLG09011531	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 60/46, Ext. diam. 65 mm, Length 165 mm	
LLG09011532	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 71/51, Ext. diam. 75 mm, Length 165 mm	1 pz.
LLG09011533	Sockets / Cones with NS joint, DURAN® tubing, Socket, NS 85/55, Ext. diam. 90 mm, Length 185 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

