

**Nalgene Oak-Ridge centrifuge tubes, PC****Thermo Scientific****Nalgeneä Oak-Ridge centrifuge tubes, PC**

Transparent. PC with PP screw cap. Can be used in suitable refrigerated or ambient centrifuges at up to 50,000xg rcf. Autoclavable. Accessory sealing cap assemblies are available for 10ml, 30ml and 50ml sizes.



Code	Description	Packaging
LLG09315691	Nalgene Oak-Ridge centrifuge tubes, PC, 10 ml, Diam. 16.1 mm, Height 81.7 mm	1 pz.
LLG09315694	Nalgene Oak-Ridge centrifuge tubes, PC, 50* ml, Diam. 28.8 mm, Height 107 mm	1 pz.
LLG09315695	Nalgene Oak-Ridge centrifuge tubes, PC, 85** ml, Diam. 38.2 mm, Height 105.7 mm	1 pz.
LLG09315692	Nalgene Oak-Ridge centrifuge tubes, PC, 28 ml, Diam. 25.4 mm, Height 101.8 mm	1 pz.
LLG09315693	Nalgene Oak-Ridge centrifuge tubes, PC, 30 ml, Diam. 25.7 mm, Height 94.5 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

