

**Sealing cap assemblies Nalgene, PP****Thermo Scientific****Sealing cap assemblies Nalgeneä, PP**

PP screw cap with Tefzel® ETFE plug liner and viton O-ring. Leak-proof at recommended maximum speed and rcf for tubes. Suitable for working with aggressive chemicals. The 38mm assembly has a PP shoulder reinforcement ring to prevent collapsing of the 250ml bottles.



Code	Description	Packaging
LLG09315734	Sealing cap assemblies Nalgene, PP, For tubes 250 ml, Screw cap diam. 38 mm	1 pz.
LLG09315732	Sealing cap assemblies Nalgene, PP, For tubes 30 ml, Screw cap diam. 20 mm	1 pz.
LLG09315733	Sealing cap assemblies Nalgene, PP, For tubes 50 ml, Screw cap diam. 24 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

