



Precision balances Type EMB

Precision balances Type EMB

The Precision balances EMB are entry-level devices in the field of laboratory balances. They can be stacked with a stack frame (optionally available) and can also be operated by battery, as an alternative to mains power, which means that they have the flexibility to be used in different locations.

- Simple and convenient 2-key operation
- Tare function facilitates formulation work
- Flat design
- Batteries included
- EMB 100-3 with ring-shaped draught shield
- Hook for underfloor weighing
- Large LCD display
- AUTO-OFF function after 3 min

Dimensions (W x D x H):

EMB 100-3/EMB1000-2:

Net weight: 170 x 240 x 39 mm

170 x 240 x 54 mm

approx. 750 g



Code	Description	Packaging
LLG06238426	Precision balances Type EMB, Type EMB 100-3, 100 g	1 pz.
LLG06235728	Precision balances Type EMB, Type EMB 1000-2, 1000 g	1 pz.
LLG09901323	Precision balances Type EMB, Type EMB 1200-1, 1200 g	1 pz.
LLG09901320	Precision balances Type EMB, Type EMB 200-2, 200 g	1 pz.
LLG09901371	Precision balances Type EMB, Type EMB 2200-0, 2200 g	1 pz.
LLG09901322	Precision balances Type EMB, Type EMB 500-1, 500 g	1 pz.
LLG09901324	Precision balances Type EMB, Type EMB 5.2K1, 5200 g	1 pz.
LLG09901372	Precision balances Type EMB, Type EMB 5.2 K 5, 5200 g	1 pz.
LLG09901321	Precision balances Type EMB, Type EMB 600-2, 600 g	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

