

**ASTM centrifuge tube, cylindrical, with conical base, borosilicate glass 3.3****ASTM centrifuge tube, cylindrical, with conical base, borosilicate glass 3.3**

Cylindrical, with conical base. Graduations and labelling in high-contrast blue enamel. Maximum loading 700 xg rcf. Internal diam. neck approx. 17 mm.

Graduated, ml: In subdivisions, ml:

0 to 0.5 0.05

0.5 to 2 0.10

2 to 3 0.20

3 to 5 0.50

5 to 10 1.00

10 to 25 5.00

25 to 100 25.0



Code	Description	Packaging
LLG09315820	ASTM centrifuge tube, cylindrical, with conical base, borosilicate glass 3.3, former standard ASTM D 96, 100 ml, Length 167 mm	1 pz.
LLG09315810	ASTM centrifuge tube, cylindrical, with conical base, borosilicate glass 3.3, acc. to ASTM D 91, 100 ml, Length 203 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

