

**Test tube shaker LLG-uniTEXER 1****LLG Labware****Test tube shaker LLG-uniTEXER 1**

Compact and **cost-efficient** test tube shaker with brushless motor drive and touch function for short-term operation. Adjustable speed control up to 4500 rpm for excellent mixing results in a very short time. Especially for small volumes.

Scope of supply: vortexer, 1 mains adapter each US, AUS, UK

Specifications

Type of movement: orbital

Orbit: 4.5 mm

Speed: 0 ... 4500 min⁻¹

Permissible ambient temperature: 5 ... 40 °C

Permissible relative moisture: 80 %

Dimensions (W x D x H): 100 x 100 x 70 mm

Weight: 0.55 kg

Power supply: 100 ... 240 V, 50/60 Hz

IP code: IP 43

Warranty: 3 years



Code	Description	Packaging
LLG06263490	uniTEXER 1 vortexer with EU-Plug	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

