



Sleeves with ribs, PTFE

Sleeves with ribs, PTFE

For use with ground glass joints. Suitable for vacuums up to 0.1mbar. Reusable and long lasting. The sealing ribs grip on the inside wall of the glass socket, sealing without the need for grease and enabling the cone to be freed easily when required.



Code	Description	Packaging
LLG09011651	Sleeves with ribs, PTFE, NS 7/16	
LLG09011652	Sleeves with ribs, PTFE, NS 10/19	
LLG09011653	Sleeves with ribs, PTFE, NS 12/21	1 pz.
LLG09011655	Sleeves with ribs, PTFE, NS 14/23	1 pz.
LLG09011656	Sleeves with ribs, PTFE, NS 19/26	1 pz.
LLG09011657	Sleeves with ribs, PTFE, NS 24/29	1 pz.
LLG09011658	Sleeves with ribs, PTFE, NS 29/32	1 pz.
LLG09011659	Sleeves with ribs, PTFE, NS 34/35	
LLG09011660	Sleeves with ribs, PTFE, NS 45/40	1 pz.
LLG09011661	Sleeves with ribs, PTFE, NS 60/46	1 pz.
LLG09011662	Sleeves with ribs, PTFE, NS 71/51	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

