


**Cryovial racks Nalgene, PC****Thermo Scientific****Cryovial racks Nalgeneä, PC**

Autoclavable. Space-saving. Choice of sizes to hold 25 or 50 cryovials of 1.2, 2.0 or 5 ml capacities. With moulded identification numbers and letters. One-handed operation possible.



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09400958 | Cryovial racks Nalgene, PC, Colour Blue, Width 197 mm, Depth 102 mm, Height 22 mm, Aperture array 5 x 5, staggered qty. | 1 pz. |
| LLG09400944 | Cryovial racks Nalgene, PC, Colour White, Width 197 mm, Depth 102 mm, Height 28 mm, Aperture array 5 x 10 qty. | 1 pz. |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

