

**Nalgene Freezing container Mr. Frosty, PC****Thermo Scientific****Nalgene Freezing container Mr. Frosty, PC**

With 120mm diameter screw cap in blue HDPE, vial holder made of white HDPE. Foam interior. Cooling rate 1°C/min. Tubes can be easily removed. Each aperture has a moulded number for identification. Stackable.



Code	Description	Packaging
LLG09400945	Nalgene Freezing container Mr. Frosty, PC, For tubes 1.0 ... 2.0 ml, Diam. 117 mm, Height 86 mm, No. of tubes 18	1 pz.
LLG09400990	Nalgene Freezing container Mr. Frosty, PC, For tubes 3.6 ml, Diam. 117 mm, Height 151 mm, No. of tubes 12	1 pz.
LLG09400991	Nalgene Freezing container Mr. Frosty, PC, For tubes 4.5 ... 5.0 ml, Diam. 117 mm, Height 151 mm, No. of tubes 12	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

