



Vertical cryobox racks Nalgene, Type 5036

Thermo Scientific

Vertical cryobox racks Nalgeneä, Type 5036

Stainless steel. With individual retainers to secure each box safely.



Code	Description	Packaging
LLG09400953	Vertical cryobox racks Nalgene, Type 5036, Type 5036 , Width 140 mm, Depth 143 mm, Height 225 mm, Height Box approx. 52 mm, Shelves 4	1 pz.
LLG09400954	Vertical cryobox racks Nalgene, Type 5036, Type 5036 , Width 140 mm, Depth 143 mm, Height 502 mm, Height Box approx. 52 mm, Shelves 9	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

