


**Floating cryovial racks Nalgene, Type 5974, PP****Thermo Scientific****Floating cryovial racks Nalgeneä, Type 5974, PP**

For 05 to 2.0ml test tubes and cryovials. Racks will float with a full load of filled tubes. With moulded identification numbers and letters.



Code	Description	Packaging
LLG09400984	Floating cryovial racks PP, Type 5974, Circular, For tubes 1.0/1.2/1.5/2.0 ml, Colour White, dim. Ø 96 mm	4 pz.
LLG09400982	Floating cryovial racks PP, Type 5974, Square, For tubes 1.0/1.2/1.5/2.0 ml, Colour White, dim. 103 x 103 x 65 mm	1 pz.
LLG09400983	Floating cryovial racks PP, Type 5974, Circular, For tubes 1.0/1.2/1.5/2.0 ml, Colour White, dim. Ø 66 mm	4 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

