

**LLG-Pressure-seal bags, PE****LLG Labware****LLG-Pressure-seal bags, PE**

Transparent. PE. With rim channel seal that can be opened and closed repeatedly.



| Code        | Description                                                            | Packaging |
|-------------|------------------------------------------------------------------------|-----------|
| LLG09404161 | Pressure-seal bags, PE, Width 100 mm, Length 150 mm, Thickness 0.05 mm | 100 pz.   |
| LLG09404165 | Pressure-seal bags, PE, Width 200 mm, Length 300 mm, Thickness 0.05 mm | 100 pz.   |
| LLG06304307 | Pressure-seal bags, PE, Width 220 mm, Length 310 mm, Thickness 0.05 mm | 100 pz.   |
| LLG09404166 | Pressure-seal bags, PE, Width 250 mm, Length 350 mm, Thickness 0.05 mm | 100 pz.   |
| LLG09404167 | Pressure-seal bags, PE, Width 300 mm, Length 400 mm, Thickness 0.05 mm | 100 pz.   |
| LLG09404160 | Pressure-seal bags, PE, Width 80 mm, Length 120 mm, Thickness 0.05 mm  | 100 pz.   |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

