



LLG-Pressure-seal bags with write on patch, PE

LLG-Pressure-seal bags with write on patch, PE

Transparent. PE. With rim channel seal that can be opened and closed repeatedly, and write-on patch.



Code	Description	Packaging
LLG09404171	Pressure-seal bags with write on patch, PE, Width 100 mm, Length 150 mm, Thickness 0.05 mm	100 pz.
LLG09404172	Pressure-seal bags with write on patch, PE, Width 120 mm, Length 170 mm, Thickness 0.05 mm	100 pz.
LLG09404177	Pressure-seal bags with write on patch, PE, Width 160 mm, Length 220 mm, Thickness 0.05 mm	100 pz.
LLG09404178	Pressure-seal bags with write on patch, PE, Width 200 mm, Length 300 mm, Thickness 0.05 mm	100 pz.
LLG09404176	Pressure-seal bags with write on patch, PE, Width 40 mm, Length 60 mm, Thickness 0.05 mm	100 pz.
LLG09404173	Pressure-seal bags with write on patch, PE, Width 70 mm, Length 100 mm, Thickness 0.05 mm	100 pz.
LLG09404179	Pressure-seal bags with write on patch, PE, Width 80 mm, Length 120 mm, Thickness 0.05 mm	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

