



Sealing tape, PTFE

Sealing tape, PTFE

For covering areas which require protection against chemicals of all types at temperatures of -150 to +260°C. Non-swelling. Easy loosening of seams even after continuous operation. Supplied in rolls as outlined below.



Code	Description	Packaging
LLG09404952	Sealing tape, PTFE, Width 12.8 mm, Length 670 cm, Wall thickness 0.1 mm	1 pz
LLG09404955	Sealing tape, PTFE, Width 24 mm, Length 670 cm, Wall thickness 0.1 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

