



## Conductivity meters Seven2Go pro S7

Mettler-Toledo Online

**Conductivity meters Seven2Go pro S7**

Portable professional model for conductivity measurement.

- New, intuitive menu structure with fast and easy T-Pad navigation
- Robust design with IP67 protection for use in harsh and demanding environments
- Intelligent Sensor Management (ISM®) functionality reduces the risk of errors
- New sunlight-proof screen guarantees screen readability even in direct sunlight
- New StatusLight gives a fast, visual confirmation showing if the meter is ready to measure or not
- Operating support for 12 languages with information shown in complete words and sentences
- Easy 1 point calibration with option of manual cell constant entry
- 3 predefined standards and 1 user defined standard
- Automatic, manual and time dependant endpoint recording
- Data memory storage for 2000 GLP data points with weatherproof Micro-USB connection for data export
- Power requirements: 4 x 1.5 V AA or 1.3 V NiMH accumulators

**Scope of supply:**

**S7-Meter:** Includes instrument, LabX direct pH, USB cable, electrode clip, base accessory, operating instructions, quick guide, test certificate, declaration of conformity, batteries

**S7-Standard kit:** As S7-Meter, but also with InLab 738-ISM and sachets of conductivity standard 1413  $\mu\text{S/cm}$  and 12.88  $\text{mS/cm}$

**S7-Field kit:** As S7-Standard kit, but also with bottles for calibration, uGoä carrying case

**S7-USP/EP kit:** As S7-Field kit, but with InLab 742-ISM instead of InLab 738-ISM

**Specifications****Conductivity****Measuring range:****Resolution:****Accuracy:**

0.010  $\mu\text{S/cm}$  ... 1000  $\text{mS/cm}$

0.001 ... 1 (auto-range)

$\pm 0.5\%$

**Temperature****Measuring range:****Resolution:**

Accuracy: -5,0 ... 105  $^{\circ}\text{C}$

0.1  $^{\circ}\text{C}$

$\pm 0.1\text{ }^{\circ}\text{C}$

TDS



**Measuring range:**

**Resolution:**

**Accuracy:**

**0.01 mg/l ... 600 g/l**

**0.01 ... 1 (auto-range)**

**±0.5 %**

**Specific resistance**

**Measuring range:**

**0.01 ... 100.0 MWcm**

**Salinity**

**Measuring range:**

**0.01 ... 42 psu**

**Conductivity Ash**

**Measuring range:**

**0.01 ... 2022 %**

**Dimensions (W x D x H): 222 x 70 x 35 mm**

Code	Description	Packaging
LLG06265256	Conductivity meters Seven2Go pro S7, Type S7-Field kit	1 pz.
LLG06265254	Conductivity meters Seven2Go pro S7, Type S7-Meter	1 pz.
LLG06265255	Conductivity meters Seven2Go pro S7, Type S7-Standard kit	1 pz.



Code	Description	Packaging
LLG06265257	Conductivity meters Seven2Go pro S7, Type S7-USP/EP kit	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

