

**Disposable spoons for foodstuffs SteriPlast, PS****Bürkle GmbH****Disposable spoons for foodstuffs SteriPlast®, PS**

Blue. The spoon for foodstuffs is available in two sizes, 2.5 ml (similar to a teaspoon) and 10 ml (similar to a tablespoon), and is ideal for sampling powders, granulates, pastes and fluids. The shade of blue that is used does not occur in a natural form in foodstuffs. Because of this, the spoons or parts of them can be quickly found and easily identified in a visual check if they are ever mixed into production. The detectable spoons have a special additive in the material that makes them visible to metal detectors or X rays. They can therefore be rejected with standard systems for checking for foreign objects, even though they are made of plastic.

- Individually packed, gamma-sterilized
- Production, assembly and packaging according to clean room class 7 (10000)
- According to EU food and FDA guidelines



Code	Description	Packaging
LLG06284560	Disposable spoons for foodstuffs SteriPlast®, PS, Visual identification, 2.5 ml pack of 100	
LLG06284561	Disposable spoons for foodstuffs SteriPlast®, PS, Visual identification, 10 ml pack of 100	100 pz
LLG06284571	Disposable spoons for foodstuffs SteriPlast, PS, Detectable , 2.5 ml pack of 100	100 pz
LLG06284572	Disposable spoons for foodstuffs SteriPlast, PS, Detectable, 10 ml pack of 100	100 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

