



## Single-use hypodermic needles, Sterican

### Single-use hypodermic needles, Sterican®

- In accordance with ISO 7864 and DIN 13097
- Thin-walled needles
- Made of stainless, chromium-nickel steel
- Smooth surface with light silicone coating
- Minimal pain upon puncture
- Transparent Luer-Lock plastic hub
- Made of Polypropylene (needle hub)
- Box of 100 pieces



Code	Description	Packaging
LLG06076115	Single-use hypodermic needles, Sterican®, Gauge 22 x 1 1/4", Diam. 0.7 mm, Colour black, Length 30 mm	100 pz
LLG06076221	Single-use hypodermic needles, Sterican®, Gauge 21 x 1 1/2", Diam. 0.8 mm, Colour green, Length 40 mm	100 pz
LLG06080870	Single-use hypodermic needles, Sterican®, Gauge 20 x 1 1/2", Diam. 0.9 mm, Colour yellow, Length 40 mm	100 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

