

**Beard Mask, White, PP, Detectable****Beard Mask, White, PP, Detectable**

This beard mask has a detectable strip on the inside. That's why it can be detected by a detector machine at any time- even if the beard mask gets into the production.

-air-permeable and breathable

-ideal for special hygienic-sensitive areas



Code	Description	Packaging
LLG06260903	Beard Mask, White, PP, Detectable, Beard mask, detectable, Colour white	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

