

**Disposable scraper for foodstuffs SteriPlast, PS****Bürkle****Disposable scraper for foodstuffs SteriPlast®, PS**

Blue. Sturdy and convenient scraper for disposable use. The scraper is suitable for scraping off, smoothing or cleaning as well as removing residue. The shade of blue that is used does not occur in a natural form in foodstuffs. Because of this, the scrapers or parts of them can be quickly found and easily identified in a visual check if they are ever mixed into production. The detectable scrapers have a special additive in the material that makes them visible to metal detectors or X rays. They can therefore be rejected with standard systems for checking for foreign objects, even though they are made of plastic.

- Individually packaged, gamma-sterilized
- Production, assembly and packaging according to Cleanroom class 7 (10000)
- According to EU food and FDA guideline



Code	Description	Packaging
LLG06284573	Disposable scraper for foodstuffs SteriPlast®, PS, Detectable, Width 80 mm, Length 200 mm	10 pz
LLG06284562	Disposable scraper for foodstuffs SteriPlast®, PS, Visual identification, Width 80 mm, Length 200 mm	10 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

