

**Moebius stirrer, stainless steel 1.4571****IKA****Moebius stirrer, stainless steel 1.4571**

Generates axial flow. For drawing the material to be mixed from above. For circulating and handling low to medium viscous liquids. Recommended for medium speed operations.



Code	Description	Packaging
LLG06265593	Moebius stirrer, stainless steel 1.4571, Type R 3000.1, Ø agitator 100 mm, Ø agitator shaft 10 mm, Length 565 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

