

**Blade stirrer, stainless steel 1.4571****IKA****Blade stirrer, stainless steel 1.4571**

Generates a radial flow. For drawing the material to be mixed from the top and bottom. For aerating of low viscosity media and for homogenizing solids in liquids. Used at medium to high speeds.



Code	Description	Packaging
LLG06265598	Blade stirrer, stainless steel 1.4571, Type R 3004, Ø agitator 30 mm, Ø agitator shaft 8 mm, Length 359 mm	1 pz
LLG06265599	Blade stirrer, stainless steel 1.4571, Type R 3004.1, Ø agitator 50 mm, Ø agitator shaft 10 mm, Length 565 mm	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

