



Disposable Syringes Injekt Solo, 2-piece

Disposable Syringes Injekt® Solo, 2-piece

- Graduation in ml for volume dosing, scale is extended over the nominal volume
- Material: Barrel: Polypropylene, Plunger: Polyethylene
- Highly transparent barrel with green plunger and black graduations for ideal readability, permanent marking
- Safe plunger backstop for easy pulling up to the maximum volume
- Silicone oil-free
- Luer Slip to put on the needle, centric or eccentric or Luer Lock conus, centric
- Latex- and PVC-free
- Sterile, single packed or non-sterile bulk packing
- Manufactured according to EN ISO 7886-1



Code	Description	Packaging
LLG06078449	Disposable syringes Injekt®Solo sterile 10 ml / Nozzle Luer, eccentric , Grad. 0,5 ml, single packed / Pack of 100	100 pz
LLG06307167	Disposable syringes Injekt®Solo sterile 2 ml / Nozzlelock, centric , Grad. 0,1 ml, single packed / Pack of 100	100 pz
LLG06800110	Disposable syringes Injekt®Solo sterile 2 ml / Nozzle Luer, centric , Grad. 0,1 ml, single packed / Pack of 100	100 pc
LLG06070497	Disposable syringes Injekt®Solo sterile 20 ml / Nozzle Luer, eccentric , Grad. 1 ml, single packed / Pack of 100	100 pc
LLG06302682	Disposable syringes Injekt®Solo sterile 5 ml / Nozzle Luer, eccentric , Grad. 0,2 ml, single packed / Pack of 100	100 pezzi



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

