

**Overhead shaker Reax 20****Heidolph****Overhead shaker Reax 20**

- Overhead Shaker for 4, 8 or 12 bottles or mixing cylinders with 160 to 270 mm height and max. 136 mm diameter
- Ideal for waste water analysis acc. to DIN 38414 Part 4
- Bottles can be easily inserted and removed by quick-release technology
- Set and continuously adjust the variable speed on the analog control knob
- Other rotation speeds available upon request, for example 0.5 to 8 rpm or 2 to 32 rpm

Specifications**Type of movement: overhead****Speed: 1 ... 16 rpm****Rotation angle: 360°****Permissible shaking weight: 30 kg****Motor rating input: 280 W****Power supply: 230 V / 50 Hz****Protection class: IP 21**

Code	Description	Packaging
LLG09839805	Overhead shakers REAX 20,for 4 bottles cap.2 ltr speed: 1 - 16 rpm	1 pz.
LLG09839806	Overhead shakers REAX 20,for 8 bottles cap.2 ltr speed: 1 - 16 rpm	1 pz.
LLG09839807	Overhead shakers REAX 20,for 12 bottles cap.2 ltr speed: 1 - 16 rpm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

