

**Filter holders, vacuum, Type 162. Stainless steel**

Sartorius Lab Instruments

Filter holders, vacuum, Type 162. Stainless steel

Specially designed for residue analysis. The stainless steel filter support produces homogeneous distribution of the residues on the membrane filter surface. Easy to use.

For 47mm or 50mm diameter filters.



Code	Description	Packaging
LLG09604202	Filter holders, vacuum, Type 162. Stainless steel, 100 ml	1 pz.
LLG09604201	Filter holders, vacuum, Type 162. Stainless steel, 500 ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

