

**Glass vacuum filter holders, glass frit**

Sartorius Lab Instruments

**Glass vacuum filter holders, glass frit**

Two compact, glass vacuum filter holders for simple particle testing and clarification of solutions. With glass frit membrane filter support and anodised aluminium clamp.

**Scope of delivery**

Size 25 mm, 30 ml: Glass vacuum filtration unit with filter support consisting of PTFE ring and glass frit.

Size 47/50 mm, 250 ml: Glass vacuum filtration unit with silicone cover and filter support consisting of PTFE ring and glass frit.



Code	Description	Packaging
LLG09604204	Glass vacuum filter holders, glass frit, Type with glass frit, For Filter size 25 mm dia., Volume 30 ml	1 pz.
LLG09604307	Glass vacuum filter holders, glass frit, Type with glass frit, For Filter size 47/50 mm dia., Volume 250 ml	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

