

**Laboratory refrigerators and freezers LR / LF series, up to +1 °C / -30 °C****Laboratory refrigerators and freezers LR / LF series, up to +1 °C / -30 °C**

- White steel exterior, stainless steel interior
- Automatic door closing
- Air forced cooling, adjustable shelves, castors, key lock on door
- Micro processor controller with digital display
- Approx. 72 hours battery back up for alarms, loggings and temperature display in case of power cut
- Prepared for GSM alarm
- Prepared for connection of 2 additional probes
- Auto cycle in case of probe failure
- Ambient temperature display
- Direct download of logged data on USB memory stick, direct upload of new software via USB memory stick
- 3-level password protected
- Battery level indication
- Display text available in different languages
- Temperature graph, integrated data logger
- RS485/RS232 interface
- Integrated memory for 30 years
- Visual and acoustic alarm
- Adjustable high/low temperature alarm
- Alarms: Power failure, probe failure, instrument failure, open door
- Contact for remote alarm, shows all alarms
- 5 year warranty on the compressor

Further models available on request.



Code	Description	Packaging
LLG06283371	Laboratory refrigerators LR 270-2, 2 x 161 l, External dim. (W x D x H) 520 x 700 x 1997 mm	
LLG06283372	Laboratory refrigerators LR 500, 515 l, External dim. (W x D x H) 620 x 860 x 1997 mm	
LLG06264456	Laboratory freezers LF 1400, 1361 l, External dim. (W x D x H) 1440 x 860 x 1997 mm	1 pz.
LLG07971187	Laboratory refrigerators LR 1400, 1361 l, External dim. (W x D x H) 1440 x 860 x 1997 mm	1 pz
LLG06283373	Laboratory refrigerators LR 660-2, 2 x 288 l, External dim. (W x D x H) 720 x 860 x 1997 mm	
LLG06283377	Laboratory freezers LF 660-2, 2 x 288 l, External dim. (W x D x H) 720 x 860 x 1997 mm	1 pz.
LLG06283375	Laboratory freezers LF 300, 346 l, External dim. (W x D x H) 520 x 700 x 1997 mm	1 pz.
LLG07971003	Laboratory refrigerators LR 300, 346 l, External dim. (W x D x H) 520 x 700 x 1997 mm	
LLG06283376	Laboratory freezers LF 500, 515 l, External dim. (W x D x H) 620 x 860 x 1997 mm	1 pz.
LLG06259589	Laboratory freezers LF 700, 618 l, External dim. (W x D x H) 720 x 860 x 1997 mm	1 pz.
LLG06283370	Laboratory refrigerators LR 100, 94 l, External dim. (W x D x H) 610 x 659 x 790 mm	
LLG06267198	Laboratory freezers LF 100, 94 l, External dim. (W x D x H) 610 x 659 x 790 mm	1 pz.
LLG07971100	Laboratory refrigerators LR 700, 618 l, External dim. (W x D x H) 720 x 860 x 1997 mm	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
 Reagents operates with
 a Certified Quality
 Management System



