

**Disk mill PULVERISETTE 13 classic line****Fritsch****Disk mill PULVERISETTE 13 *classic line***

Ideal for intermittent and continuous fine grinding within the medium particle size range of soft to hard, tough and temperature-sensitive solids in the areas of mining and metallurgy, glass industry, ceramics industry, soil research, rocks and soils.

Max. feed size: 20 mm

Throughput: up to 150 kg/h

Final fineness: 0.1 to 12 mm

Supply requirements: 400 V/3~, 50 Hz, 1790 W

Dimensions (W x D x H): 440 x 870 x 400 mm

Net weight: 140 kg



Code	Description	Packaging
LLG09739001	Disc mill "pulverisette 13" without milling disc 400V/3-, 50-60Hz, 2050W	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

