



## Reagent solutions for photometers Aqualytic

### Reagent solutions for photometers Lovibond®

Solution with 15 ml or refill bottle with 100 ml.



Code	Description	Packaging
LLG06058068	Reagent solutions for photometers Lovibond®, For Chlorine, DPD No. 1, Colour blue	1 pz.
LLG07510822	Reagent solutions for photometers Lovibond®, For Chlorine, DPD No. 1*, Colour blue	1 pz.
LLG06057943	Reagent solutions for photometers Lovibond®, For Chlorine, DPD No. 1, Colour green	1 pz.
LLG06057944	Reagent solutions for photometers Lovibond®, For Chlorine, DPD No. 3, Colour red	1 pz.
LLG06203285	Reagent solutions for photometers Lovibond®, For Chlorine, DPD No. 3*, Colour red	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

