

**pH electrode A 7780****SI Analytics****pH electrode A 7780**

Glass shaft, Silamid®-reference system, with 1 m or 1.5 m (A 7780 IDS) fixed cable.

Specifications

Temperature range: -5 ... 80 °C

pH range: 0 ... 14 pH

Zero point: pH=7.0 ±0.3

Taper length / Dia.: 120mm / 12mm

Electrolyte: Gel

Diaphragm: 3 x ceramic

Membrane shape: Sphere

Membrane resistance at 25 °C: 250 MOhm



Code	Description	Packaging
LLG06265619	pH electrode A 7780, Type A 7780 IDS, Temperature sensor NTC 30 kW, Connector Digital	1 pz.
LLG06235995	pH electrode A 7780, Type A 7780 1M-DIN-ID, Temperature sensor Pt 1000, Connector DIN + 4-mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

