



Refractive Index standards

Reagecon

Refractive Index standards

Certified standards with values expressed in Refractive Index (R.I) units nD.

- ISO 17025 accredited with an uncertainty of measurement up to ± 0.00014 nD units at 20 °C
- 12 month shelf life
- Aqueous based
- Compatible with any brand of refractometer
- Convenient high quality 15 ml HDPE dropper bottle



Code	Description	Packaging
LLG06265877	Refractive Index standards, Refractive Index at 20 °C 1.34325nD	1 pz.
LLG07658877	Refractive Index standards, Refractive Index at 20 °C 1.34782nD	1 pz.
LLG06265878	Refractive Index standards, Refractive Index at 20 °C 1.37233nD	1 pz.
LLG06258375	Refractive Index standards, Refractive Index at 20 °C 1.38115nD	1 pz.
LLG06265879	Refractive Index standards, Refractive Index at 20 °C 1.39032nD	1 pz.
LLG06265142	Refractive Index standards, Refractive Index at 20 °C 1.40978nD	1 pz.
LLG06258376	Refractive Index standards, Refractive Index at 20 °C 1.42009nD	1 pz.
LLG06265880	Refractive Index standards, Refractive Index at 20 °C 1.44193nD	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

