



## UV analysis lamp HP-UVIS

### UV analysis lamp HP-UVIS®

For UV analysis without a darkroom. Inclusive UV shield. High-pressure mercury lamp emits very intense radiation at 366nm for fluorescence evaluation. Specially selected 254nm filters enable optimal contrast. Minimal footprint in spite of simultaneous observation points for two 200mm x 200mm plates. Inclined plate angle gives comfortable viewing.

Dimensions (W x D x H) 325mm x 290mm x 480mm. Weight: 11.5kg. Supply requirements 230V.



Code	Description	Packaging
LLG09539360	UV analysis lamp HP-UVIS®, Type HP-UVIS®	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

