


**Test tubes, Durham, soda lime glass****Mariefeld****Test tubes, Durham, soda lime glass**

Durham test tubes are mainly used in microbiology to detect gas produced by micro organisms.

These special tubes are placed upside down in bigger test tubes.

- Made of soda lime glass
- With straight rim and with round bottom
- Ungraduated



Code	Description	Packaging
LLG06266183	Test tubes, Durham, soda lime glass, Diam. 6 mm, Wall thickness 0.7 mm, Height 35 mm	250 pz.
LLG06266184	Test tubes, Durham, soda lime glass, Diam. 6 mm, Wall thickness 0.7 mm, Height 50 mm	350 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

