

**Dissolved oxygen meter FiveGo F4**

Mettler-Toledo Online

Dissolved oxygen meter FiveGo F4

Portable dissolved oxygen meter for field work.

- Measurement made easy: A large, well-structured display, intuitive button arrangement and simple menu ensure measurements can be performed in just a few clicks
- Ergonomic design for one handed operation: Thanks to its ergonomic design, handling the FiveGo F4 is a comfortable experience. The foldable stand provides optimal viewing on flat surfaces

Scope of supply:**F4-Meter:** Meter only**F4-Standard kit:** Meter with LE621 sensor**F4-Field kit:** Meter with LE621 sensor, carrying case**Specifications****Measuring range (O₂ Saturation):****Resolution:****Accuracy:** 0.1 ... 400 %

0.1 %

±1 %

Measuring range (Temperature):**Resolution:****Accuracy:** 0 ... 50 °C

0.1 °C

±0.3 °C

Measuring range (O₂ concentration):**Resolution:****Accuracy:** 0.01 ... 45.00 mg/l (ppm)

0.01 mg/l

±1 %

Temperature compensation: ATC and MTC**Calibration:** 2 points, 100 % and 0 %**Memory size:** 200 measurements, current calibration**Dimensions (W x D x H):** 188 x 77 x 33 mm**IP code:** IP 67 (meter + sensor)

Code	Description	Packaging
LLG06280120	Dissolved oxygen meter FiveGo F4, Type F4-Field kit	1 pz.
LLG06280118	Dissolved oxygen meter FiveGo F4, Type F4-Meter	1 pz.
LLG06280119	Dissolved oxygen meter FiveGo F4, Type F4-Standard kit	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

