

**Steel Bottle Cart KKR-832/01 and KKR 820/1****Gorr****Steel Bottle Cart KKR-832/01 and KKR 820/1**

Steel bottle cart with solid rubber tyres Ø 200 mm and powder coated black-silver.

Bottle holder with chain protection. Wheels: metal rim, roller bearing.



Code	Description	Packaging
LLG09148304	Steel Bottle Cart KKR-832/01 and KKR 820/1, For 10 l, Type KKR-820/01 for 1 flask	1 pz.
LLG09148303	Steel Bottle Cart KKR-832/01 and KKR 820/1, For 33 l, Type KKR-832/01 for 1 flask	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

