



Accessories for Knife Mill GM 200/300

Accessories for Knife Mill GM 200/300



Code	Description	Packaging
LLG09738506	Accessories for Knife Mill GM 200/300, Type GM 300, Gravity lid, plastic, autoclavable and scratch-resistant	1 pz.
LLG09738507	Accessories for Knife Mill GM 200/300, Type GM 300, Gravity lid, plastic, with overflow channels, autoclavable and scratch-resistant	1 pz.
LLG09739210	Accessories for Knife Mill GM 200/300, Type GM 200, Gravity lid, PP for grinding container of autoclavable plastic and PP	1 pz.
LLG09739211	Accessories for Knife Mill GM 200/300, Type GM 200, Gravity lid, PP for grinding container of glass or steel	1 pz.
LLG09739218	Accessories for Knife Mill GM 200/300, Type GM 200, Gravity lid, PP with overflow channels for grinding container of autoclavable plastic and PP	1 pz.
LLG09739209	Accessories for Knife Mill GM 200/300, Type GM 200, Gravity lid, PP with overflow channels for grinding container of glass or steel	1 pz.
LLG09738510	Accessories for Knife Mill GM 200/300, Type GM 200, Grinding container, 1 litre, plastic, autoclavable	1 pz.
LLG09738501	Accessories for Knife Mill GM 200/300, Type GM 300, Grinding container, plastic, transparent, autoclavable and scratch-resistant, 5L, including knife holder	1 pz.
LLG09738502	Accessories for Knife Mill GM 200/300, Type GM 300, Grinding container, plastic, transparent, autoclavable and scratch-resistant, 5L, without knife holder	1 pz.
LLG09739213	Accessories for Knife Mill GM 200/300, Type GM 200, Grinding container, PP, 1L	1 pz.
LLG09739248	Accessories for Knife Mill GM 200/300, Type GM 200, Grinding container, stainless steel, 1L	1 pz.
LLG09738503	Accessories for Knife Mill GM 200/300, Type GM 300, Grinding container, stainless steel, autoclavable, 5L, including knife holder	1 pz.
LLG09738504	Accessories for Knife Mill GM 200/300, Type GM 300, Grinding container, stainless steel, autoclavable, 5L, without knife holder	1 pz.
LLG09739254	Accessories for Knife Mill GM 200/300, Type GM 200, Knife, pure titanium for heavy-metal free grinding	1 pz.
LLG06267312	Accessories for Knife Mill GM 200/300, Type GM 300, Knife, stainless steel (full metal blade and cylinder), autoclavable	



Code	Description	Packaging
LLG06267311	Accessories for Knife Mill GM 200/300, Type GM 200, Knife, stainless steel with serrated blades	1 pz
LLG06267313	Accessories for Knife Mill GM 200/300, Type GM 300, Knife, stainless steel with serrated blades, autoclavable	1 pz
LLG06267314	Accessories for Knife Mill GM 200/300, Type GM 300, Knife, titanium-niob coating, autoclavable, for grinding without heavy metal contamination	
LLG09739222	Accessories for Knife Mill GM 200/300, Type GM 200, Reduction lid, PP for volume reduction to 0.5L for grinding container of autoclavable plastic or PP	1 pz.
LLG09738509	Accessories for Knife Mill GM 200/300, Type GM 300, Scraper	1 pz.
LLG09739212	Accessories for Knife Mill GM 200/300, Type GM 200, Scraper	1 pz.
LLG07930431	Accessories for Knife Mill GM 200/300, Type GM 200, Stainless steel knife for GM 200	1 pz.
LLG09738508	Accessories for Knife Mill GM 200/300, Type GM 300, Standard knife, stainless steel, autoclavable	1 pz.
LLG09738505	Accessories for Knife Mill GM 200/300, Type GM 300, Standard lid, plastic, transparent, autoclavable and scratch-resistant	1 pz.
LLG07930430	Accessories for Knife Mill GM 200/300, Type GM 200, Standard lid, PP for grinding container of autoclavable plastic, PP, glass or steel	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

